

SCOUT DS®

Noninvasive Diabetes Screening for Type 2 and Prediabetes


NO BLOOD • NO FASTING • NO WAITING



Urgency of early diabetes detection.

Diabetes is a costly global epidemic due in part to the micro- and macrovascular complications that begin years before diagnosis.^{1,2} Additionally, much of diabetes can be prevented or delayed if detected early.³ However, despite comprehensive screening guidelines, diagnosis of diabetes typically doesn't occur until 5–9 years post onset when 50% of patients have one or more complications.^{4,5}

NATURAL HISTORY OF TYPE 2 DIABETES^{1,2}




Better metrics for early diabetes detection.

Increases in glucose dysregulation and oxidative stress cause measurable changes in early diabetes-related biomarkers. The SCOUT Diabetes Score is derived via measurements of fluorescence and reflectance from these biomarkers found in the skin of the forearm, including:

- Measures of glycemc exposure and oxidative stress including Advanced Glycation End Products (AGEs) and connective tissue cross-links
- Measures of metabolic change including NADH and FAD
- Measures of microvascular changes

A proprietary algorithm transforms measurements from these biomarkers, corrects for skin tone, and reports a SCOUT Diabetes Score.



SCOUT DS® Noninvasive diabetes screening.

Benefits

Noninvasive

No blood draw or biohazards

No fasting

Allows opportunistic testing any time

Fast, immediate results

Facilitates patient counseling at point of service

Easy to use

Three-step procedure with factory calibration and on-board controls

Superior sensitivity^{6,7}

Contributes to earlier disease detection

No test strips to handle

Just enter pre-paid test key



SIMPLE TEST PROCEDURE

- 1 Enter patient age and gender
- 2 Place arm in cradle
- 3 Results in as little as 80 seconds

SCOUT DS® shows superior detection

of abnormal glucose tolerance versus either FPG or A1C...

TEST ⁶	THRESHOLD	SENSITIVITY (%) [95% CI]	FALSE POSITIVE RATE (%) [95% CI]	ABSOLUTE SENSITIVITY DIFFERENCE	RELATIVE SENSITIVITY DIFFERENCE
SCOUT	SCOUT DIABETES SCORE = 50	75.2 [70.8 – 79.2]	42.1 [37.4 – 47.0]	–	–
FPG	ADA 5.6 MMOL/L (100 MG/DL)	56.0 [51.1 – 60.7]	23.7 [19.9 – 28.1]	19.2%	34%
	WHO 6.1 MMOL/L (110 MG/DL)	37.6 [33.1 – 42.4]	7.4 [5.2 – 10.3]	37.6%	100%
A1C	ADA 5.7%	72.5 [67.9 – 76.6]	43.5 [38.8 – 48.3]	2.7%	4%
	IEC 6.0%	46.8 [42.0 – 51.6]	15.1 [11.9 – 18.8]	28.4%	61%

... and of dysglycemia (HbA1c \geq 6.0%) versus Random Capillary Glucose...


TEST ⁷	THRESHOLD	SENSITIVITY (%) [95% CI]	FALSE POSITIVE RATE (%) [95% CI]	ABSOLUTE SENSITIVITY DIFFERENCE	RELATIVE SENSITIVITY DIFFERENCE
SCOUT	SCOUT DIABETES SCORE = 50	83.8 [76.7 – 89.1]	32.8 [27.4 – 38.7]	–	–
RCG	6.1 MMOL/L (110 MG/DL)	33.1 [25.7 – 41.4]	19.1 [14.8 – 24.3]	50.7%	254%

... and of diabetes (HbA1c \geq 6.5%) versus Random Capillary Glucose.


TEST ⁷	THRESHOLD	SENSITIVITY (%) [95% CI]	FALSE POSITIVE RATE (%) [95% CI]	ABSOLUTE SENSITIVITY DIFFERENCE	RELATIVE SENSITIVITY DIFFERENCE
SCOUT	SCOUT DIABETES SCORE = 50	98.0 [89.3 – 99.6]	43.6 [38.5 – 48.8]	–	–
RCG	6.1 MMOL/L (110 MG/DL)	44.9 [31.9 – 58.7]	20.9 [17.0 – 25.5]	53.1%	218%

Interpretation of results.

SCOUT DS®



The SCOUT Diabetes Score reference method is the 2 hour 75g Oral Glucose Tolerance Test (OGTT).



* Based on a prospective study of 471 subjects at risk for type 2 diabetes by American Diabetes Association criteria and 267 subjects with previously diagnosed type 2 diabetes.

Who should be screened?

In general, patients aged 40 or older should be tested approximately every three years.

If any of the listed risk factors apply, they should be tested earlier and/or more often.

- Overweight
- Physical inactivity
- Insulin resistance
- Family history of diabetes
- High risk ethnic group
- Vascular disease
- Hypertension
- Dyslipidemia
- Complications associated with diabetes
- Gestational diabetes
- Delivery of a macrosomic infant
- Polycystic ovary syndrome

1. Ramio-Halsted B.A., et al. *Primary Care*. 1999;26:771–789.
2. Nathan D.M. *New England Journal of Medicine*. 2002;347:1342–1349.
3. DREAM (Diabetes Reduction Assessment with ramipril and rosiglitazone Medication) Gerstein H.C. et al. *The Lancet*. 2001;344:1343–1350.
4. Harris M.I., et al. *Diabetes Metab Res Rev*. 2001;16:230–236.
5. Manley S.M., et al. *Diabetes Res*. 1990;13:1–11.
6. Sensitivity, specificity, and associated 95% confidence intervals for SCOUT DS, the FPG test, and the A1C test from the ENGINE trial. Results based on an analysis set of 471 subjects, all at risk for type 2 diabetes by American Diabetes Association criteria.
7. Sensitivity, specificity, and associated 95% confidence intervals for SCOUT DS and RCG test from the GREECE trial. Results based on an analysis set of 409 subjects, all at risk for type 2 diabetes by American Diabetes Association criteria.



800 Bradbury Drive SE, Suite 217
Albuquerque, NM 87106, USA
+1 (505) 272-7023

www.veralight.com